

EXHIBIT C

From Page No. 2

Purpose - To label synthetic probe ChtA176 with ³²P by kinasing for Chlamydia assay development.

Reagents -

Synthetic probe, ChtA176 mm 320-39 1.28 OD/ml 64 µg/ml
 10X Kinase buffer, from EB + mm, aliquot obtained (391:18)
 T4 Kinase, BRL, 10 u/µl Lot 62111
 γ-³²P-ATP 3000 Ci/mmol #G551, 10 µCi/µl
 Teflon, acid-base treated, from Mark H.
 5 m NaCl DEB NB 157:48
 1% SDS diluted from 20% (DK)
 1 M Tris, pH 8.2 (DK)
 100 mM EE L/N 500:38
 4 M LiCl from DK
 Proteinase K 20 mg/ml 333:11
 Glycogen HL-1 ~40 mg/ml in 10% EtOH (DK)
 Phenol from DK ()
 Chloroform
 10% TCA
 BSA, 10 mg/ml from BRL Lot# 40416

Procedure

1. Follow same procedure used (391:18-20)
 Kinase reaction included 1 µl of ChtA176 (64 ng)

Results

Column rinse in 0.3 M NaCl and H₂O:

$$\approx (9.6 \times 10^4) \left(\frac{3300}{5} \right) = 6.3 \times 10^7 \text{ cpm}$$

Straight count = $(6.2 \times 10^4) \left(\frac{200}{1} \right) \left(\frac{80}{10} \right) = 6.2 \times 10^7 \text{ cpm}$ > why so low?

Input dpm = $2.2 \times 10^8 \text{ dpm}$

EtOH super = $(1.1 \times 10^3) \left(\frac{2400}{5} \right) = 5.3 \times 10^5 \text{ cpm}$

Witnessed & Understood by me,

Joann KOP

Date

Invented by

Mary E. Harper

Recorded by

Date

To Page No. 26

TITLE ^{32}P -Kinasing of ChtA176

Pr j ct N _____

Bo k No. 391

26

From Page No 25

Results, cont

$$\text{TGA pptable cpm} = (1.2 \times 10^4) \left(\frac{200}{1} \right) \left(\frac{50}{10} \right) = 1.2 \times 10^7 \text{ cpm in } 64 \text{ ng}$$

$$\text{specific act.} = 1.9 \times 10^8 \text{ cpm}/\mu\text{g}$$

$$\text{Final cpm} \quad (1.8 \times 10^5) \left(\frac{100}{2} \right) = 9 \times 10^6 \text{ cpm} \quad 75\% \text{ recovery}$$

$$\text{in } 100 \mu\text{l H}_2\text{O} \quad 9 \times 10^4 \text{ cpm}/\mu\text{l}$$

Diluted total sample to 300 μl 0.02% SDS $\rightarrow 3 \times 10^4 \text{ cpm}/\mu\text{l}$ final
store at 4°C

To Page No. _____

Witnessed & Understood by me,

Joann Kop

Date _____

Invented by _____

Date _____

Recorded by _____

Mary E. Harper